

MEMO

TO: Clients

FROM: CCLS – Specimen Referral Center

DATE: October 4, 2022

SUBJECT: Fentanyl with Metabolite Confirmation, Random, Urine

Effective October 25, Mayo will be making changes to the following test:

FENTANYL WITH METABOLITE CONFIRMATION, RANDOM, URINE, MAYO CODE FENTU (LAB1186)

Explanation: A method enhancement will result in the following test changes.

Current Specimen Required
<p>Supplies: Urine Tubes, 10 mL (T068) Collection Container/Tube: Clean, plastic urine collection container Submission Container/Tube: Plastic, 10-mL urine tube Specimen Volume: 10 mL Specimen Minimum Volume: 2.1 mL Collection Instructions: 1. Collect a random urine specimen. 2. No preservative. Additional Information: 1. No specimen substitutions. 2. No STATS are accepted for this procedure. 3. Submitting less than 10 mL will compromise our ability to perform all necessary testing.</p>

New Specimen Required
<p>Supplies: Sarstedt Aliquot Tube, 5 mL (T914) Collection Container/Tube: Clean, plastic urine collection container Submission Container/Tube: Plastic, 5-mL tube Specimen Volume: 5 mL Specimen Minimum Volume: 1 mL Collection Instructions: 1. Collect a random urine specimen. 2. No preservative. Additional Information: 1. No specimen substitutions. 2. No STATS are accepted for this procedure. 3. Submitting less than 5 mL may compromise our ability to perform all necessary testing.</p>

Current Specimen Stability		
Specimen	Temperature	Time
Urine	Refrigerated (preferred)	14 days
	Ambient	72 hours
	Frozen	14 days

New Specimen Stability		
Specimen	Temperature	Time
Urine	Refrigerated (preferred)	28 days
	Ambient	21 days
	Frozen	28 days

Please feel free to contact the Specimen Referrals Center with any questions at 320-251-2700 Ext. 57320