

MEMO

TO: Clients

FROM: CCLS – Specimen Referral Center

DATE: 8/24/2021

SUBJECT: Inhibin A

Effective September 21st, 2021, Mayo will be making changes to the following test:

INHIBIN A, TUMOR MARKER, SERUM, MAYO CODE INHA (LAB1791)

The preferred specimen type and specimen stability are changing to the following:

| Current Specimen Required | |
|----------------------------|-----------|
| Collection Container/Tube: | |
| Preferred: | Red Top |
| Acceptable: | Serum Gel |

| New Specimen Required | |
|----------------------------|-----------|
| Collection Container/Tube: | |
| Preferred: | Serum Gel |
| Acceptable: | Red Top |

| Current Specimen Stability | | |
|----------------------------|--------------------------|---------|
| Specimen | Temperature | Time |
| Serum | Refrigerated (preferred) | 7 days |
| | Frozen | 90 days |

| New Specimen Stability | | |
|------------------------|--------------------------|---------|
| Specimen | Temperature | Time |
| Serum | Refrigerated (preferred) | 14 days |
| | Frozen | 90 days |

Please feel free to contact the Specimen Referrals Center with any questions at 320-251-2700 Ext. 57320